


Calcium Microscope



Nikon Ti-E 3-Camera Calcium Widefield Microscope

Calcium Mic - STEP BY STEP INSTRUCTIONS CONTENTS PAGE

Page 5	-----	STEP 1	System On
Page 10	-----	STEP 2	Lens And Focus
Page 15	-----	STEP 3	Change From Eyes To Camera View
Page 17	-----	STEP 4	Setting Up Initial Live View
Page 25	-----	STEP 5	Optimising Your Camera Settings (When You Don't Have Enough Signal)
Page 36	-----	STEP 5	Optimising Your Camera Settings (When You Have Too Much Signal)
Page 40	-----	STEP 5	Optimising Your Camera Settings (When There Is Too Much Background Noise – Use Averaging)
Page 43	-----	STEP 5	Optimising Your Camera Settings (Now Do Step 5 For Every Channel You Want To Use)
Page 44	-----	STEP 6	Acquisition Settings (Save To File)
Page 47	-----	STEP 6	Acquisition Settings (Order Of Acquisition Tabs)
Page 49	-----	STEP 6	Acquisition Settings (Lambda Λ Tab – Channel Selection)
Page 52	-----	STEP 6	Acquisition Settings (Z-stack)
Page 57	-----	STEP 6	Acquisition Settings (Large Image)
Page 60	-----	STEP 6	Acquisition Settings (XY Positions)
Page 62	-----	STEP 6	Acquisition Settings (Time)
Page 66	-----	STEP 7	At The End Of Your Session (Save And Shut Down Procedures)
Page 71	-----	MORE	More Advanced Instructions (Reuse Pervious Camera Settings)
Page 73	-----	MORE	More Advanced Instructions (Time Measurement)
Page 75	-----	!!!	KEY INSTRUCTIONS TO REMEMBER AND FOLLOW


Before Using The Facility...

To use the facility the individual must undergo training...

- Tour of the facilities
- T & C agreements
- First training session (2 hours)
- Second training session with user's samples (2 hours)
- Additional training sessions may be required
- Additional lens installation training (at a later date)
- Access to booking system: <http://ppms.eu/kcl-wohl>

STEP BY STEP INSTRUCTIONS

The rest of this document will take you from focusing on your sample to optimising your image to what to do at the end of your session.



Additional
Information
Slides

The blue slides contains additional information you might find helpful.

STEP BY STEP INSTRUCTIONS

STEP 1

System ON

The Equipment

DAPI/GFP/RFP/Far-red LAMP


Brightfield

DG4 FURA-2 LAMP

Camera 7C

Camera 7B

Camera 7A



System On

1. Switch on the Calcium Microscope depending on the function you need **(consult the flowchart on the wall)** and follow the numbered switches.
2. **Make sure the stage is empty before turning on the Microscope switch at the right, far back.**
3. **ALWAYS** login to NIS-Elements Software before loading any sample, this checks if all systems are connected.





Software ON : Choose A Camera

NIS elements will ask you to choose a camera to use:

ANDOR Neo/Zyla
Switch number 7A
(7A and 7B if you are doing two-camera imaging) Is for Brightfield / DIC / DAPI / GFP / RFP / Far-Red / CFP / YFP / FURA-2

Nikon DS-Fi3
Switch number 7C
Is for coloured imaging.





Software ON : Initial Settings

1. After loading NIS elements, please check you have the standard panels loaded (circled in red).
2. Top Left, in OC Panel, in Eyes, select Brightfield (eyes).

OC Panel

Eyes

- DAPI (eyes) | GFP (eyes) | mCherry (eyes)
- CFP (eyes) | **Brightfield (eyes)**

Camera

- DAPI | GFP | DsRed
- FarRed | Brightfield | ColourCameraBF

DIC Camera

- DIC | DIC shade correction

Fura-2 Camera

- 340 Fura 2 | 380 Fura 2 | 470nm green
- 405nm green

Two Camera

- GFP_mCherry | GFP_mCherry_59022

OC Panel Camera

- YFP | CFP | CFP FRET

ND Acquisition

Save to File | Record Data...

Custom Metadata

Order of Experiment | Timing...

Time | XY | Large Image | λ | z

Setup

Opt. Conf.	Name	Comp. Color	Focus Of...
<input checked="" type="checkbox"/> FarRed	FarRed		0
<input type="checkbox"/> DsRed	DsRed		0
<input checked="" type="checkbox"/> GFP	GFP		0
<input type="checkbox"/> DAPI	DAPI		0
<input type="checkbox"/> DIC shade correction	DICsc		0

Close Active Shutter during Filter Change

Use Ratio | Define Ratio...

Advanced >>

Load | Save | Remove | 1 time loop | Run now

Neo Settings

Format: No binning

Auto Exposure: 100 ms

1 frame for Fast Timelapse

Readout Mode: Rolling shutter

Readout Rate: 200 MHz

Dynamic Range: 12-bit & Gain 1

Sensor Mode: Overlap

Limit Maximum FPS to: 25 < 10.0 FPS

Spurious Noise Filter

Temperature: -31.0 °C

Commands

F3000 Pad | Filters, Shutters and Switchers

F3000

LUTs


No Image

STEP BY STEP INSTRUCTIONS

STEP 2

Lens And Focus

Checking Lens for Damage and Cleanliness



Make sure you are on brightfield eyes and white light is coming through. Check the lens for any damage by putting the toggle beneath the eyepiece to the **LEFT** and zoom on the lens using the button here.

ALWAYS put the toggle back to the **RIGHT**
after you're done and before you start to focus.

Cleaning OIL Lens

1. Wipe away any excess oil with dry lens tissue (normally after your session)
 2. Wrap lens tissue around your finger and soak up some ethanol and clean lens from centre outwards (REPEAT 3 TIMES)
 3. Clean once more with dry lens tissue
- During your session if you are switching between oil/dry lenses, just wipe away excess oil on the lens and your slide, lower your objectives down, before switching.

Changing Lens

ALWAYS

Lower the lenses as far as they can go before inserting the stage and/or clicking on another lens.

Nosepiece



To lower lens...

- Rotate wheel clock-wise
- Bottom = this number no longer reduces
- Be careful - zero doesn't mean it's at the bottom

Lower the lenses between changing slides.



Focus



ALWAYS put the toggle
back to the **RIGHT**
after you're done and
before you start to focus.

DRY LENS ONLY (20X Air and 40X Air)

- Use Wheel and joystick on the control pad to adjust and focus.
- Make sure you are on the specimen.
- Focus indicator comes when it detects a reflective surface. (Do not rely on numbers as the lowest setting is not always set to ZERO, to set the stage XZ to 0 μ m press here.)
- When the green light comes on STOP moving up immediately (lower it back down to when the light just flash on if needed, otherwise the focal plane might be passed).
- Switch coarse XZ to fine (Z button), coarse XY to fine (twist the joystick)
- Select one of the "EYE" options in the software.
- Find focus using eye piece.



OIL LENS

- **Focus indicator does NOT work for oil/water lens.**
- Wheel upwards on coarse ONLY until lens comes in contact with oil.
- Switch coarse XZ to fine, coarse XY to fine (twist the joystick)
- Find focus using eye piece.
- Find focus on camera.


STEP BY STEP INSTRUCTIONS

STEP 3

Change from Eyes to Camera View


From this point forward we will be adjusting things on the computer screen so we need to work in camera mode instead of looking down the eye piece.

Remain on STOP setting.



Click on a 'camera' option


This moves from eyes to camera.



STEP BY STEP INSTRUCTIONS

STEP 4


Setting Up Initial Live View



STEP 1: Good initial settings to view samples in live is 100ms exposure and gain 1, this can be adjusted later depending on how fluorescent your sample is.


Step 2: In camera mode, live screen comes up automatically when you click this.

Now you're in live, to visualise your image you need to adjust your contrast.



To adjust contrast in LUT Panel

1. Click on auto contrast
2. Starting to see something on screen, but it might be out of focus.
3. Hover over your live image and adjust the mouse wheel slowly to change focus.
4. Slide this bar to the top and adjust contrast line further if needed.



TI Pad x Neo Settings x

Format: No Binning

Auto Exposure: 5 ms

1 frame for Fast Timelapse

Readout Mode: Rolling shutter

Readout Rate: 560 MHz

Dynamic Range: 12-bit & Gain 1

Sensor Mode: Overlap

Limit Maximum FPS to: 25 < 100.2 FPS

Spurious Noise Filter

Temperature: -30.3 °C

Commands

F3000 Pad x Filters, Shutters and Switchers x

0 100 11 [%]

F3000

LUTs x

ON V W X Y Z

Auto Scale

94 G: 1.00 110

0 1000 2000 3000 4000

OC Panel x

Eyes

DAPI (eyes) GFP (eyes) mCherry (eyes)

CFP (eyes) Brightfield (eyes)

Camera

DAPI GFP DsRed

FarRed Brightfield ColourCameraBF

DIC Camera

DIC DIC shade correction

Fura-2 Camera

340 Fura 2 380 Fura 2 470nm green

405nm green

Two Camera

GFP_mCherry GFP_mCherry_59022

CFP YFP Camera

YFP CFP CFP FRET

ND Acquisition x ND Sequence Acquisition x

Experiment: ND Acquisition

Save to File Record Data...

Custom Metadata

Order of Experiment Timing...

Time XY Large Image

Setup

Opt. Conf.	Name	Comp. Color	Focus Of...
FarRed	FarRed		0
DsRed	DsRed		0
GFP	GFP		0
DAPI	DAPI		
DIC shade correction	DICsc	Brightfield	0

Close Active Shutter during Filter Change

Use Ratio Define Ratio...

Advanced >>

Load Save Remove 1 time loop Run now

LUTs

Additional Information Slides...

LUTs And Contrast

LUTs explained.

- Changing LUTs to visualise your sample better does not change your raw data (signal intensity).
- This means when you open your image again in e.g. ImageJ, LUTs will not be applied.
- If you change LUTs before quantitative analysis, it is recommended to save the LUTs and apply it to all comparable images.

Why do we need to change the contrast?

The camera in this microscope captures shades from 0 to 4000, initially the LIVE window shows you all these different shades, but the signal from this sample only reach roughly 2500 therefore we only need to work within the 0 to 2500 range.

Auto contrast brings the contrast into a range for you to better visualise your sample.

1 Click on this

Initial live window

Auto contrast applied

LUTs x Auto Scale G: 1.00 0 4095

LUTs x Auto Scale G: 1.00 0 2824

Auto Contrast

Reset LUTs

LUTs can be saved in a separate file and applied to other images.

Continuous auto contrast (there may be flickering while in live view)

Auto Contrast

Reset LUTs

LUTs in more detail...

Continuous auto-contrast while in live view – can cause a flickering effect.

Auto-contrast


Delete all contrast adjustments

To keep analysis consistent, you can copy and paste LUTs across different captured images.

Makes dim targets brighter while keeping bright target the same.

Controls the Y axis log graph

Controls the X axis contrast, this doesn't change the captured raw data.



The spinning disk does not over saturate, but on other microscopes, this is for viewing saturation, complementary colour recommended. Once an area is saturated, it loses any intensity value information.

The camera is a black and white camera. You can assign any colour combination to your captured image.

Drag and change the threshold of this histogram for better visualisation.

X axis - relative to #-bit camera (e.g. 16-bit CCD camera gives 65,536 different intensity values, of which you should not go over 50,000). Y axis – log intensity scale

Fit the histogram to this space

STEP BY STEP INSTRUCTIONS

STEP 5

Optimising your camera settings

...when you don't have enough signal

Changing LUTs such as contrast does not change your raw data (i.e. it does not affect how much signal your sample is giving off)

To optimise your signal, we need to optimise the camera settings.


How much signal is enough signal?

That really depends on what you want to measure...

Hover your mouse over the background to measure background intensity read it here

Then hover over something you want to measure, cell body for example and the intensity read out will change accordingly..

If the difference in intensity between background and target is under 50, its probably not good enough for accurate analysis, if you are counting cells by eye, ~200 difference is recommended, if you want to do any automated thresholding or counting, over 800 difference is recommended.




FarRed 0.32 $\mu\text{m}/\text{px}$ Mono 12bit: 2160 x 2560 pixels [1287, 1369] Mono: 100

To increase signal intensity...


Step 1: Reduce exposure to ~10ms


Step 2: Change Gain1 to Gain4 (this amplifies your signal drastically)



Gain 1

Gain 4







Gain (everything in this box) need to be the same across ALL the channels you are going to use.

I recommend using your weakest channel to set your gain, then for other channels if you have too much signal, just reduce your exposure time.

Gain 1





Gain 4





Remember to click on auto-contrast

Gain 4

Gain 4





After auto-contrast





Now check background intensity and target intensity again, difference of ~ 200 is good for counting by eye but not for automated counting. Adjust the exposure to optimise further.



Increase exposure time to 20ms...



Auto-contrast and check intensity difference again.



>1000 difference between background and target intensities, usable for auto thresholding.



STEP BY STEP INSTRUCTIONS


STEP 5

Optimising your camera settings

...when you have too much signal


How to tell if you have too much signal?


Oversaturated pixel



When your LUTs graph is filled up like this with a peak at the end, you may be oversaturating your sample.

Turn on your oversaturating indicator by selecting complementary colour in the drop down.





Oversaturation!

Oversaturation means the camera is picking up too much signal, and the camera can no longer determine the actual intensity of your signal, it just knows that sample is 'bright'.

This can cause you problems during analysis, because you won't have intensity information.

To fix oversaturation...

1. Reduce exposure time.
2. You can also reduce Gain - but if you can't change this because you are limited by another channel (see rule on page 32) then just reduce exposure time further.


STEP BY STEP INSTRUCTIONS

STEP 5

Optimising your camera settings

...when there is too much background NOISE


use AVERAGING.




background
NOISE

Averaging

- Takes multiple images and averages them out.
- Beware 8X averaging increase your imaging time by 8 folds.




Step 1: Zoom In and then go on LIVE to better see if your image is noisy and need averaging.



Step 2: Determine if you need averaging and how many time you need to average your image.

Happy with camera settings... remember to **SAVE!!!**

- ! means there's been a change in your camera settings, such as gain and averaging, exposure time is automatically updated and saved.



SAVE current camera settings for Far-Red channel into the Far-Red optical configuration by clicking on the  make sure you click on the right button.

STEP BY STEP INSTRUCTIONS

STEP 5

Optimising your camera settings

...now do STEP 5 for every channel you want to use.

STEP BY STEP INSTRUCTIONS


STEP 6

Acquisition Settings

... Save to File

Set Up File Path (SAVE)





- 1) Check Save to file option
- 2) Go to Browse and select DATADRIVE1:
- 3) Create/find your folder, set up new folder for this session if needed.
- 4) Recommended file name:
Experiment_Name_Date_001
- 5) Finished, anytime you press "Run now" a new file will be automatically saved in the nd2 format.



Run VS Capture

Automatically saved if
'Save to File' is ticked

Not automatically saved



Box will appear
and you can save
or discard image.

After 'Run now' if you make changes to your image (such as adding ROIs) then save this 'new' image by going to file and 'SAVE AS' so you don't overwrite your raw data.


STEP BY STEP INSTRUCTIONS

STEP 6

Acquisition Settings

... Order or acquisition tabs

Fastest acquisition tab sequence



The Software will prioritise the tab on the RIGHT.

Do not put 'Large Image' tab on the right.


Tick the box for all acquisition functions you want to use.

STEP BY STEP INSTRUCTIONS


STEP 6

Acquisition Settings


... Lambda (channel selection)



Opt	Name	C	Focus Of...
<input checked="" type="checkbox"/>	FarRed		
<input type="checkbox"/>	DsRed		0
<input type="checkbox"/>	GFP		0
<input type="checkbox"/>	DAPI		
<input type="checkbox"/>	DIC correction		Brightfield 0



Add and select all the channels you want to use



NIS-Elements AR [Current user: Chen] - [Cells_001.nd2*]

File Edit Acquire Calibration Image ROI Binary Measure Reference Macro View Devices Window Applications JOBS Help

OC Panel x

Eyes: DAPI (eyes), GFP (eyes), mCherry (eyes), CFP (eyes), Brightfield (eyes)

Camera: DAPI, GFP, DsRed, FarRed, Brightfield, ColourCameraBF

DIC Camera: DIC, DIC shade correction

Fura-2 Camera: 340 Fura 2, 380 Fura 2, 470nm green, 405nm green

Two Camera: GFP_mCherry, GFP_mCherry_59022

CFP YFP Camera: YFP, CFP, CFP FRET

ND Acquisition x ND Sequence Acquisition x

Experiment: ND Acquisition

λ: []

Save to File

Path: E:\Chen [Browse...]

Filename: Cells_001.nd2 [Record Data...]

Custom Metadata

Order of Experiment Timing...

Time [] XY [] Large Image [] λ [] z []

Setup

Opt. Conf.	Name	Comp. Color	Focus Of...
<input checked="" type="checkbox"/>	FarRed	[Magenta]	0
<input type="checkbox"/>	DsRed	[Red]	0
<input type="checkbox"/>	GFP	[Green]	0
<input type="checkbox"/>	DAPI	[Blue]	0
<input type="checkbox"/>	DIC shade correction	Brightfield	0

Close Active Shutter during Filter Change

Use Ratio [Define Ratio...]

Advanced >>

Load Save Remove 1 time loop Run now

Cells_001.nd2*

ND Progress

Experiment overall progress:

λ: []

Time elapsed: 0:00:00 Time remaining: 0:00:00

Experiment Status:

- Lambda

Detail Info

Events...

Pause Refocus Finish Abort

F3000 Pad x Neo Settings x

Format: No Binning

Auto Exposure: 20 ms

1 frame for Fast Timelapse

Readout Mode: Rolling shutter

Readout Rate: 560 MHz

Dynamic Range: 12-bit & Gain 4

Sensor Mode: Overlap

Limit Maximum FPS to 25 < 50.0 FPS

Spurious Noise Filter

Temperature: -30.3 °C

Commands

F3000 Pad x Filters, Shutters and Switchers x

0 100 11 [%]

F3000

LUTs x

G: 1.00 4095

0 1000 2000 3000 4000

FarRed 0.32 μm/px Mono 12bit: 2160 x 2560 pixels [1062, 1797] Mono: 0

FarRed 0.32 μm/px Mono 12bit: 2160 x 2560 pixels [N/A]

Full Screen Measurement Basic widefield* Neo (90 deg) ActivateDocument("Cells_001.nd2"); Plan Apo VC 20x DIC N2 (0.32 μm/px @ 2560 x 2160) XY=[-20.960, 1.124]mm, Z=[] Logged in as chenliang (48m)

Click on 'Run now' to start imaging.

STEP BY STEP INSTRUCTIONS

STEP 6

Acquisition Settings

... Z stack



Z stack Basic Options


Set top and bottom: use mouse wheel to focus and define the exact range of your Z stack.

Set Middle: use mouse wheel to find the mid point of your Z focus and set equal distance above and below the focal plane. (Useful if your sample is symmetrical along the Z axis)

Asymmetrical: find focal plane and then set different distances above and below. (Useful for  like cells)

Step Size

- You can set step size or number of steps.
- Use recommended step size to capture all the information.
- Fewer steps than the recommended is called under-sampling and you may lose information.
- More steps than the recommended is called over-sampling (typed in here) which is required for 3D deconvolution.



The screenshot shows the Z stack software interface. At the top, there are checkboxes for 'Time', 'XY', 'Large Image', and 'Z'. Below these are three icons for setting the Z stack: 'Set Top and Bottom', 'Set Middle', and 'Asymmetrical'. A 3D model of a blue rectangular prism represents the Z stack, with labels for 'Top' (1959.05 abs), 'Bottom' (1948.15 abs), and 'Range' (10.90 μm). Below the 3D model are input fields for 'Step' (0.9 μm), 'Bottom' (1948.15 μm), and 'Top' (1959.05 μm). To the right, there are fields for 'Range' (10.90 μm), 'Relative Positions' (Top: -4.17 μm, Bottom: -15.08 μm), and 'Direction' (Bottom to Top selected). At the bottom, there are buttons for 'Load', 'Save', 'Remove', '1 time loop', and 'Run now'.


Choose an option to set your Z stack.



1

Tick to select

2



N/A

N/A

4

In Live mode use mouse wheel to define top/bottom/focus of your sample.

Watch the numbers to orientate if you're going up or down.

Reset

3

Reset

5 Set step size...Click (for recommended step size)

Step: 0.9

μm

0.9 μm

14

Steps

Range:

Bottom: N/A

Or type in here to Over/Under-sample

N/A μm

Relative Positions:

Top: N/A μm

Bottom: N/A μm

Close Active Shutter during Z Movement

Direction: Bottom to Top

Top to Bottom

Advanced >>

Load

Save

Remove

1 time loop





Run now

6


Run now to image your Z stack

Z stack View Modes





In 3D rendering you have different rendering modes.



STEP BY STEP INSTRUCTIONS

STEP 6

Acquisition Settings

... Large Image



Time XY Large Image λ Z

1

Tick to select, Lambda must be ticked too if not using Z stack.

Set scan area

3


X axis x Y axis fields
5.0 x 9.0 mm
 Pattern

2

In Live mode, move to the centre of your field of view.

4

Stitch Use overlapping image edge to stitch
 Progressive Registration



5

Overlap: 15 % 15% overlap minimum!

Close Active Shutter during Stage Movement

Run now to image your Large Image

6


Load Save Remove Advance Run Z Corr 1 time loop Run now

Click on this to fit your whole large image into the window

Do not put Large Image tab on the far right.

Y=1

X=3




STEP BY STEP INSTRUCTIONS

STEP 6

Acquisition Settings

... XY Positions



ALWAYS DELETE ALL PREVIOUS POSITIONS BEFORE YOU START!

ALWAYS REMOVE ALL WHEN YOU'VE FINISHED AND UNCHECK XY TAB

Time XY Large Image

Points Move Stage to Selected Point + Add

Point Name	X [mm]	Y [mm]	Z [μm]	PFS
<input checked="" type="checkbox"/> #1	-20.960	1.124	1962.78	N/A
<input checked="" type="checkbox"/> #2	-21.826	2.523	1961...	N/A
<input type="checkbox"/>				

Include Z Relative XY

Close Active Shutter during Stage Movement

1

6

3

4

5

2

Add that field of view to your positions.

In Live, drag and move to another field of view. (repeat 2-4 until you have all your positions)

Run now : the microscope will scan and move to the next position as fast as it can.

In Live, focus on your field of view.

STEP BY STEP INSTRUCTIONS

STEP 6

Acquisition Settings

... Time



Time Schedule

Time XY Large Image λ Z

2

Phase	Interval	Duration	Loops
<input checked="" type="checkbox"/> #1	1 sec	10 sec	11
<input type="checkbox"/>		msec	
		sec	
		min	
		hour(s)	

Add a time-lapse
Interval (how long between each scan)
Duration (how long the overall time-lapse)
Loops (automatically calculated for you)

Close Active Shutter when idle Perform Time Measurement (0 ROIs)
 Switch Transmitted Illuminator off when Idle

Events... Advanced >>

Load Save Remove 1 time loop Run now

1

In Live, focus on your field of view.

3

Run now

NIS-Elements AR [Current user: Chen] - [Cells_005.nd2 T(1/11): 731.260 msec FPS overall: 1.01]

File Edit Acquire Calibration Image ROI Binary Measure Reference Macro View Devices Window Applications JOBS Help

Search

OC Panel x

Eyes

- DAPI (eyes) GFP (eyes) mCherry (eyes)
- CFP (eyes) Brightfield (eyes)

Camera

- DAPI GFP DeRed
- FarRed Brightfield ColourCameraBF

DIC Camera

- DIC DIC shade correction

Fura-2 Camera

- 340 Fura 2 380 Fura 2 470nm green
- 405nm green

Two Camera

- GFP_mCherry GFP_mCherry_59022

CFP YFP Camera

- YFP CFP CFP FRET

ND Acquisition x ND Sequence Acquisition x

Experiment: ND Acquisition

T: A:

Save to File

Path: E:\Chen Browse...

Filename: Cells_006.nd2 Record Data...

Custom Metadata

Order of Experiment Timing...

Time XY Large Image A Z

Time Schedule

Phase	Interval	Duration	Loops
<input checked="" type="checkbox"/> #1	1 sec	10 sec	11
<input type="checkbox"/>			

Close Active Shutter when idle Perform Time Measurement (0 ROIs)

Switch Transmitted Illuminator off when Idle

Events... Advanced >>

Load Save Remove 1 time loop Run now

TI Pad x Neo Settings x

Format No Binning

Auto Exposure 20 ms

1 frame for Fast Timelapse

Readout Mode Rolling shutter

Readout Rate 560 MHz

Dynamic Range 12-bit & Gain 4

Sensor Mode Overlap

Limit Maximum FPS to 25 < 50.0 FPS

Spurious Noise Filter

Temperature: -30.3 °C

Commands

F3000 Pad x Filters, Shutters and Switchers x

0 100 11 [%]

F3000

LUTs x

88 G: 1.00 2413

1000 2000 3000 4000

5

5 6 7 8 9 10 11

T>


FarRed 0.32 µm/px Mono 12bit: 2160 x 2560 pixels [N/A]

Plan Apo VC 20x DIC N2 (0.32 µm/px @ 2560 x 2160)

XY=[-21.454, 2.914]mm, Z=[] Logged in as chenliang (1h3m)

Full Screen Measurement Basic widefield Neo (90 deg) LUTs_AutoScale0; Play Sequence

You can play your
time-lapse loop and
make a video.



Continue with your imaging...

STEP BY STEP INSTRUCTIONS

STEP 7

At the end of your session

... Shut down procedure

STEP 1 : Clean Up and Close Software

The screenshot shows the NIS-Elements AR software interface. The main window displays a brightfield image of cells. The left sidebar contains various toolbars and panels, including 'Eyes', 'Camera', and 'ND Acquisition'. The bottom status bar shows technical details like 'Basic widefield', 'Neo (90 deg)', and 'Plan Apo VC 20x DIC N2 (0.32 μm/px @ 2560 x 2160)'.

1 STOP

2 Put it in Brightfield Eyes

3 LOWER the objectives
REMOVE your sample
CLEAN the lens

4 Close Software


4 Close Software

Your data should already be automatically saved, check your DATADRIVE1 folder if you want to be sure before shutting off the software.

STEP 2 : Transfer Data to WCIC Shared Network Drive

1

Open file
Find and cut/copy your
saved data in
DATADRIVE1




2

Right click to open
WCIC Shared Drive
in new window


3

You need to login to this pop-up window, user name is
normally: kclad\k number
DO NOT click on remember my credentials



4

In the Network drive, open the
Nikon_TiE_Twocamera folder



5

Paste data into your
personal folder

- **DO NOT USE USBs ON ANY COMPUTERS IN THE MICROSCOPE ROOMS!**
- The Shared Drive can be accessed from the workstations (where you can use USBs), or your personal computers, from there please BACK-UP your data.

If you have trouble connecting to the shared drive, or need to re-map the network drive, please follow the instructions in this link or contact us.

https://f59fa4a0-6222-454a-9101-2ba099e49b58.filesusr.com/ugd/0c5e54_9b32bb45a8fc450bb3e65e9c79de3b44.pdf

KING'S
College
LONDON

Shared Network Drive

Wohl Cellular Imaging Centre

Slide 1 – Normal login
Slide 2 – When you can't see the WCIC shared drive icon
Slide 3 – WCIC connected but you can't access your folder


Normally to connect to shared drive, just follow Step 1 and 2.

Step 1 – In File Explorer click on wcic

Step 2 – login to shared drive

If you cannot see 'WCIC' the last user may have disconnected it.

STEP 3 : Check the booking schedule!



If someone is booked on within 2 hours, leave the system on.

If no one is using the system within

2 hours, shut down the system in

REVERSE numerical order.

If the incubator is on... shut down as usual but do not turn the incubator off, another user might be pre-warming the system.

STEP BY STEP INSTRUCTIONS

More advanced instructions

...reuse previous camera settings



What if you are imaging similar samples and want to **re-use camera settings** you've optimised before...

- 1) In NIS Elements software, open a previous image with camera settings you want to mimic.
- 2) Right click on the image once it's open
- 3) Select reuse camera settings
- 4) Be aware, this uploads camera settings only, acquisition setting (Z-stack, large image etc.) will not be reloaded



Time measurement

Time measurement (found under the measurement tab) can be used after time-lapse imaging to determine the rate of sample bleaching or recovery.

Removing Background (Changes to ROIs can be implemented at workstations/FIJI after your session)



Use the ROI options on the right of your captured image. This option allow you to subtract any background auto-fluorescence.



After any ROI changes, the new image will not be automatically saved. Form a good habit of "Save As" the new image with a new file name reflecting what changes you made to the raw data.

Just remember...

Make sure the stage is **EMPTY** before turning on the microscope switch at the right, far back.

Always put the toggle back to the **RIGHT** after you're done and before you start to focus.

Always **LOWER** the lenses as far as they can go before inserting the stage and/or clicking on another lens.

Always switch from **COARSE** to **FINE** when oil / water is in contact with lens.

GAIN need to remain the **SAME** for all Channels!!!

Acquisition panel **XY** positions: always **DELETE ALL** previous positions before you start !

DO NOT USE USBs on any computers in the microscope rooms!

If you need any help, please contact:

George Chennell (07771926760)

or

Chen Liang (07883166321) via WhatsApp

Happy Imaging!